



Be Right™



Digital titrator cartridge, sodium hydroxide, 8.00 N

Товар #:

1438101

Опасно

RUB Цена:

Связаться с Хах

Доступно



For determination of acid concentration by Digital Titrator titration. Digits x 0.01 = meq

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Спецификации

Диапазон измерений:	10 - 4000 meq/L Acid
Количество тестов:	40 - 100
Наименование метода:	Phenolphthalein End Point
Номер метода:	8233 (Base)
Параметр:	Acid-Base
Платформа:	Digital Titrator
Титрант:	8 N NaOH