



Be Right™



Reagent set, manganese, HR, periodate oxidation, 25 mL

Товар #:

2243200

Опасно

RUB Цена:

Связаться с Хах

Доступно



For determination of Manganese by the Periodate Oxidation Method, Hach Method 8034. EPA Approved for wastewater reporting (digestion required) - Federal Register, 44, (116) 34193 (June 14, 1979). Range: 0.1 - 20.0 mg/l Mn. Sample Size: 25 mL. Set Contains: Buffer Powder Pillows (Citrate type) and Sodium Periodate Powder Pillows for 100 tests.

Спецификации

Диапазон измерений:	0.1 - 20.0 мг/л Mn
Инструмент:	DR/3000, DR/2000
Количество тестов:	100
Наименование метода:	Periodate Oxidation
Номер метода:	8034
Платформа:	Powder Pillows
Соответствует требованиям EPA:	да