



Be Right™



Acid reagent powder pillows, 25 mL, pk/1000

Товар #: 104228
RUB Цена: Свяжитесь с Хах
Позвоните, чтобы узнать дату доставки

Опасно



For determination of high range silica by the silicomolybdate method (Hach method 8185) Range: 1-100.0 mg/L SiO₂. Also used in the iodometric method (Hach method 10100) for the determination of high range hypchlorite (50-150 g/L as Cl₂). Pack of 1000 powder pillows.

Спецификации

Наименование метода: Silicomolybdate
Платформа: Powder Pillows