

# Reagent, hydrazine, AccuVac® ampules, pk/25

Товар #: 2524025 Опасно  
RUB Цена: Связаться с Хах  
Доступно

For determination of Hydrazine by the p-Dimethylaminobenzaldehyde method, using AccuVac® Ampules. Hach Method 8141. Range: 4 - 600 µg/L. Pack of 25 ampules.

## Fast

AccuVac® Ampules provide quick results. Just snap, fill and read - saving you valuable time.

## Easy

No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

## Accurate

Optical quality evacuated vials assure precise measurements and reagent stability.

## Use them anywhere you need results.

Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

---

## Спецификации

|                      |   |
|----------------------|---|
| Диапазон измерений:  | 4 - 600 мкг/л N <sub>2</sub> H <sub>4</sub> |
| Инструмент:          | All except DR/850,<br>DR/820, PC II         |
| Количество тестов:   | 25  |
| Наименование метода: | p-Dimethylaminobenzaldehyde                 |
| Номер метода:        | 8141  |
| Параметр:            | Hydrazine                                   |
| Платформа:           | Accuvac                                     |