



Be Right™



Calibration kit, Stablcal turbidity standards, 2100AN / AN IS turbidimeter, 500 mL bottles

Товар #: 2659500
RUB Цена: Свяжитесь с Хах
Позвоните, чтобы узнать дату доставки

Опасно



Stablcal Standards Calibration Kit, for 2100AN and 2100AN IS Laboratory Turbidimeters. USEPA accepted for calibration purposes. Stablcal Stabilized Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. Kit includes 500 mL bottles of <0.1, 20, 200, 1000, 4000 NTU standard and a sealed vial of 7500 NTU standard.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Спецификации

Концентрация: < 0.1, 20, 200, 1000, 4000, 7500 NTU
Объем / Размер упаковки: set of 6