



Be Right™



## Reagent powder pillows, silver 1, pk/50

Товар #:

2293566

Опасно

RUB Цена:

Связаться с Хах

Доступно



Reagent for determination of Silver by the Cadion 2B Colorimetric Method. Hach Method 8120. Range: 0.02 - 0.70 mg/L Ag. Sample Size: 50 mL. Pack of 50 Pillows. Also available in reagent set 2296600.

### Спецификации

Количество тестов:

50

Наименование метода:

Colorimetric

Платформа:

Powder Pillows