



Be Right™



Reagent, DPD total chlorine, test 'n tube vials, pk/50

Товар #:

2105645

Опасно

RUB Цена:
Доступно

Связаться с Хах



For total chlorine determination by the DPD Method. Test 'N Tube Reagent Vials contain pre-measured DPD reagent. Range: 0.09 - 5.00 mg/L Cl₂. Sample Size: 10 mL

Pack of 50 vials.

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Boric Acid free!

Спецификации

Диапазон измерений:	0.09 - 5.00 мг/л Cl ₂
Инструмент:	All except PC II
Количество тестов:	50
Наименование метода:	DPD
Номер метода:	10101
Параметр:	Chlorine, Total
Платформа:	Powder Pillows