



Be Right™



Chlorine Reagent Set (Total), Ultra Low Range DPD

Товар #:

2563000

Опасно

RUB Цена:

Связаться с Хах

Доступно



For determination of ultralow Total Chlorine concentrations by the DPD Method. Hach Method 8370. USEPA accepted for drinking water and wastewater analysis. Range: 2 - 500 µg/L Cl₂. Reagents are ampuled for stability. Reagent set includes ULR Chlorine Buffer Solution (1.5 mL ampules), DPD Indicator Solution for ULR Chlorine (1.5 mL ampules) and Blanking Reagent for 20 tests.

Boric Acid free!

Спецификации

Диапазон измерений:	2 - 500 мкг/л Cl ₂
Инструмент:	All except DR 900, PC II
Количество тестов:	20
Наименование метода:	DPD
Номер метода:	8370
Параметр:	Chlorine
Соответствует требованиям EPA:	да