



Be Right™



## Reagent set, hardness (total), HR, digital titrator

Товар #:

2448100

Опасно

RUB Цена:

Связаться с Хах

Доступно



For the determination of Total Hardness, 100-4000 mg/l CaCO<sub>3</sub>, EDTA Method.

Reagent set for the determination of Total Hardness by EDTA Digital Titrator titration (Hach method 8213). Range: 100-4,000 mg/L as CaCO<sub>3</sub>. Set contains: Buffer Solution (Hardness 1), ManVer® 2 Hardness Indicator Powder Pillows and 0.800 M EDTA Digital Titrator Cartridge for 40 - 100 tests.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Спецификации

Диапазон измерений:	100 - 4000 мг/л CaCO <sub>3</sub>
Количество тестов:	100
Наименование метода:	EDTA Titration
Номер метода:	8213
Платформа:	Digital Titrator