

Buffer solution pillows, respirometric BOD nutrient, pk/50

Товар #: 2962266 Опасно
RUB Цена: Свяжитесь с Хах
Доступно

For use with BODTrak II Apparatus. Meets Standard Methods 5210 D requirements. Use one pillow per sample bottle. Included with BODTrak II Apparatus (2952400) and BODTrak II Reagent Set with Bottles (2952440).

Спецификации

Formula:	Nutrients & Buffer
Package settings:	Powder Pillows
Области применения:	BODTrak II
Объем / Размер упаковки:	50 pcs
	1 Pillow/Bottle
Описание:	Nutrients & Buffer