

# Standard solution, potassium, 100 mg/L as K (NIST), 500 mL

Товар #: 2351749  
RUB Цена: Связаться с Хах  
Доступно

## Ion Selective Electrode (ISE) standard solutions for laboratory analyses

For accuracy checks in potassium determinations

### Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

## Спецификации

Package settings:	Bottle
Концентрация:	100 &nbsp;мг/л К
Объем / Размер упаковки:	500 mL
Описание:	Potassium Standard Solution
Параметр:	Potassium
Платформа:	Раствор
Специальные указания:	Observed shelf life after opening, following recommended storage and use: 3 months
Срок службы:	2 years
Условия хранения:	10 - 25 °C (protect from light)