



Be Right™



## Standard solution, sulfite, 15 mg/L as SO<sub>3</sub> (equivalent), 500 mL

Товар #: 2408449  
RUB Цена: Свяжитесь с Хах  
Доступно

Sulfite equivalent standard solution. Prepared with stabilized Sodium Thiosulfate. Equivalent to 15 mg/L as SO<sub>3</sub>. 500 mL Bottle

---

### Спецификации

Концентрация: 15 mg/L  
Описание: Sulfite Equivalent Solution as SO<sub>3</sub><sup>2-</sup> (NIST)  
Платформа: Раствор