



Be Right™



## Calibration set, Stablcal turbidity standards, 2100P portable turbidimeter, 100 mL bottles

Товар #:

2659410

Опасно

RUB Цена:

Связаться с Хах

Доступно



Four ready to use, stabilised Formazine Primary standard solutions in bottles. Kit includes 100 mL of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU .

### Спецификации

Концентрация:

< 0.1, 20, 100, 800 NTU

Объем / Размер упаковки:

set of 4