



Be Right™



## SINGLET Single-Use pH Standards Combo, pH 7.00 & 10.01 (10 each)

Товар #: 2769820  
RUB Цена: Свяжитесь с Хах  
Доступно

### Singlet™ Solution Packs

Combi pack with Buffer Solutions pH 7.00 and pH 10.01, ready for use, closed airtight in separate sealed pouches. Accuracy  $\pm 0.02$  at 25°C.

For pH meter standardization. Sealed foil pouches of pH standardization buffers. pH is  $\pm 0.02$  at 25°C

#### Ideal for on-site analyses

Sealed pouches are easy to transport - no chance of spillage.

#### Easy to Use

Just open the pouch and insert the electrode.

#### Fresh and reliable

Foil pouches prevent contamination and air exposure.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

---

### Спецификации

Color coded?:	да
Package settings:	Singlet
pH:	7.00, 10.01
Концентрация:	7,00 pH и 10,01 pH
Объем / Размер упаковки:	pk
Описание:	Buffer Solution, pk/20 total (10 each)
Прослеживаемость по NIST:	да
Специальные указания:	Discard solution after use
Срок службы:	24 months
Условия хранения:	10 - 25 °C (protect from light)