

Буферный раствор pH 10.00, Сертификат анализа (загружается с сайта), 4 x 250 мл

Товар #: LZW9472.99
RUB Цена: Свяжитесь с Хах
Доступно

pH Buffer standard solutions for laboratory and field analyses

Буферный раствор pH 10.00, Сертификат анализа (загружается с сайта), 4 x 250 mL.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Спецификации

Color coded?:	нет
Package settings:	Bottle
pH:	10.01 ±0.02 pH at 25°C
Region:	Global
Объем / Размер упаковки:	1 L (4 x 250 mL)
Описание:	pH Buffer Calibration reagent
Прослеживаемость по NIST:	да
Специальные указания:	Observed shelf-life after opening, following recommended storage and use: 3 months
Срок службы:	12 months
Условия хранения:	10 - 25 °C (protect from light)