



Be Right™



Reagent set, phosphorus digestion

Товар #:

2459000

RUB Цена:

Связаться с Хах

Доступно

Опасно



For acid persulfate digestion of samples for Total Phosphorus (Phosphate) determination by Hach photometric methods. Requires additional reagents for measurement of the reactive phosphorus generated by the digestion. Set contains Potassium Persulfate Powder Pillows, Sodium Hydroxide Solution (5.0 N), and Sulfuric Acid (5.25 N). Approximately 100 tests.

Спецификации

Инструмент:	All
Количество тестов:	100
Номер метода:	8190
Требуется разложение:	Yes