



**Be Right™**



## Filter holder, easy pressure syringe, PALL

**Товар #:** 246800  
**RUB Цена:** Свяжитесь с Хах  
**Доступно**

Accepts 25 mm filter discs for fast filtration of 1-10 mL of sample.

**Fits most standard syringes with Luer-type fittings**

**Made of Delrin plastic with type-316 stainless steel membrane support screen**

**Autoclavable**

**Seal plugs included**

**Single filter PN 246801.**

---

### Спецификации

Внешний диаметр:	25 mm
Единицы измерения:	1 each

---

### Требуемые аксессуары

- Filter, membrane, pore size 0.45  $\mu$ m, diameter 25 mm, 100/pk (Item 2514101)