



Iron color disc test kit, model IR-18A

Товар #:

146500

Опасно

RUB Цена:

Связаться с Хак

Доступно



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Range:	0.02 - 1.00 mg/L Fe
Smallest increment:	0.02
Вес с упаковкой:	2
Диапазон измерений:	NA
Источник питания:	NA
Количество тестов:	100
Модель:	IR-18A
Наименование метода:	Color disc/ 1,10 Phenanthroline
Образец тестового набора:	IR-18A
Параметр:	Iron, medium range - As Fe**
Платформа:	Color Disc
Тип корпуса:	B