



Be Right™



## Cyanuric acid test kit, model CY-3

Товар #: 185102  
RUB Цена: Свяжитесь с Хах  
Поставляется в течение 2 недель

Опасно



A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Спецификации

Smallest increment:	20
Вес с упаковкой:	2
Диапазон измерений:	20 - 100 &nbsp;&nbsp;&nbsp;мг/л
Количество тестов:	100
Наименование метода:	Turbidimetry
Области применения:	Pool and Spa
Образец тестового набора:	CY-3
Параметр:	Cyanuric Acid
Тип корпуса:	D