



Iron and manganese color disc test kit, model IR-20

Товар #:

146300

Опасно

RUB Цена:

Связаться с Хак

Доступно



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3. Smallest I

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Диапазон измерений: 0 - 4,0 мг/л Fe

Наименьшие шаги приращения: 0.2

Области применения: Water Quality

Платформа: Color Disc