



Be Right™



Test kit, hardness & iron color disc, model HA-77

Товар #: 202300
RUB Цена: Свяжитесь с Хах
Позвоните, чтобы узнать дату доставки

Опасно



Hardness & Iron Water Quality Test Kit is the kit to use if you need better accuracy than that provided by a test cube. Utilizing drop count titration and color disc colorimetry, this kit contains 100 tests.

The kit to use for more resolution in a visual test

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Области применения: Water Quality

Платформа: Color Disc