



Be Right™



Test kit, corrosion management, model CM-1

Товар #: 2495700
RUB Цена: Свяжитесь с Хах
Позвоните, чтобы узнать дату доставки

To reduce lead and copper levels in distribution, service and domestic water lines, suppliers must stabilize aggressive water by adjusting pH, temperature, and hardness. Analysts can verify the success of their stabilization by screening samples with this rugged, portable test kit.

Uses drop count titration and color disc methods

Maintain corrosion control of aggressive distribution system waters to meet Lead and Copper Rule Requirements

Establish baseline conditions to monitor treatment or disinfection strategy changes

Calculate Langelier and Agressive Indexes

Спецификации

Количество тестов:	100
Области применения:	Corrosion Management
Платформа:	Color Disc