



Be Right™



Test kit, iron & pH color disc , model IR-23

Товар #: 256200
RUB Цена: Свяжитесь с Хах
Доступно

Utilizing the simple, accurate color disc method, this Iron & pH Water Quality Test Kit provides accurate, on-site measurement of two important parameters to help maintain water quality within the desired range. Range: 0-10 mg/L/ 4-10 pH unit.

Uses color disc to provide accurate, on-site measurement of two important parameters to help maintain water quality within the desired range

Спецификации

Области применения: Water Quality
Платформа: Color Disc