



Be Right™



## Test kit, boiler treatment control, model BTC-2, drop count titration

Товар #: 2350400  
RUB Цена: Свяжитесь с Хах  
Позвоните, чтобы узнать дату доставки

Опасно



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. This kit provides testing capability for sulfite, phosphate, and pH, offering an economical means of monitoring common boiler water treatment chemicals. Sulfite is generally added to boiler water to reduce dissolved oxygen levels, while phosphate additions help reduce scale formation. pH levels are monitored to ensure corrosion control.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus to perform analysis in the field**

---

### Спецификации

Диапазон измерений:	7.4 - 9.6 pH Units
Области применения:	Boiler Water
Платформа:	Color Disc