



Be Right™



Test kit, molybdate, model MO-LR

Товар #: 2359300
RUB Цена: Свяжитесь с Хах
Доступно

Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Спецификации

Range:	0 - 3.0 mg/L Mo
Smallest increment:	0.25
Вес с упаковкой:	2
Диапазон измерений:	NA
Источник питания:	NA
Количество тестов:	100
Модель:	MO-LR
Наименование метода:	Color disc/ Ternary complex
Образец тестового набора:	MO-LR
Параметр:	Molybdate - As Molybdenum
Платформа:	Color Disc
Тип корпуса:	D