



Be Right™



Test kit, color, model CO-1

Товар #: 223400
RUB Цена: Свяжитесь с Хах
Доступно

The Hach color test reads the apparent color created by dissolved substances that can indicate industrial, agricultural, or natural pollution.

Color disc kit reads apparent color created by dissolved substances that can indicate industrial, agricultural, or natural pollution

Dual range measurement

Спецификации

Range:	0 - 100 units
Smallest increment:	5 units
	25
Вес с упаковкой:	1.0
Диапазон измерений:	0 - 100 units
Источник питания:	NA
Количество тестов:	No Limit
Модель:	CO-1
Наименование метода:	Color disc/ APHA Platinum – Cobalt Standard
Образец тестового набора:	CO-1
Параметр:	Color - As color units
Платформа:	Color Disc
Примечание:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Тип корпуса:	D