



Be Right™



Test kit, manganese color disc, model MN-PAN

Товар #: 2350800
RUB Цена: Связаться с Хак
Поставляется в течение 2 недель

Опасно



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Range:	0 - 0.70 mg/L Mn
Smallest increment:	0.05
Вес с упаковкой:	3
Диапазон измерений:	NA
Источник питания:	NA
Количество тестов:	50
Модель:	MN-PAN
Наименование метода:	Color disc/PAN
Образец тестового набора:	MN-PAN
Параметр:	Manganese - As Mn
Платформа:	Color Disc
Тип корпуса:	D