



Be Right™



Manganese test kit, model MN-5

Товар #: 146700
RUB Цена: Свяжитесь с Хах
Доступно

Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Smallest increment:	0.1
Вес с упаковкой:	2
Диапазон измерений:	0 - 3.0 mg/L Mn
Количество тестов:	100
Наименование метода:	Color disc/Cold periodate
Образец тестового набора:	MN-5
Параметр:	Manganese - As Mn
Платформа:	Color Disc
Тип корпуса:	D