



Be Right™



Test kit, glutaraldehyde, model GT-1

Товар #: 2587200
RUB Цена: Свяжитесь с Хах
Доступно

From oil fields to cooling towers, hospitals to food plants, poultry farms to dental offices, the Glutaraldehyde Test Kit can help monitor disinfection processes. The test kit measures from 0.5 to 4000 ppm glutaraldehyde.

Helps monitor disinfection processes

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Smallest increment:	0.5
Вес с упаковкой:	8.3
Диапазон измерений:	1 - 2800 mg/L multiple ranges
Количество тестов:	100
Наименование метода:	Color disc
Образец тестового набора:	GT-1
Параметр:	Glutaraldehyde
Платформа:	Color Disc
Тип корпуса:	Q