



Be Right™



Test kit, RapidSilver™ visual

Товар #: 2674500
RUB Цена: Свяжитесь с Хах
Позвоните, чтобы узнать дату доставки

Test Kit Silver, Model RapidSilver™. Range 0-50 ppb. Silver As Ag Visual Comparison.

Provides fast, accurate information on silver discharges

Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss

Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples

Спецификации

Smallest increment:	5 мкг/л
Вес с упаковкой:	4
Диапазон измерений:	5 - 50 ppb Ag
Количество тестов:	100
Наименование метода:	Visual Comparison
Образец тестового набора:	RapidSilver™
Параметр:	Silver-As Ag
Тип корпуса:	F