

## Test kit, chloride, model CD-DT

Товар #: 2063500  
RUB Цена: Связаться с Хах  
Позвоните, чтобы узнать дату доставки

The single-parameter Model CD-DT uses the Hach Digital Titrator with Mercuric Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-8,000 mg/L. Contains approximately 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

---

### Спецификации

Smallest increment:	0.1-20
Вес с упаковкой:	4
Диапазон измерений:	10 - 8,000 mg/L Cl <sup>-</sup>
Количество тестов:	100
Наименование метода:	Digital Titrator/ Mercuric nitrate
Образец тестового набора:	CD-DT
Платформа:	Digital Titrator
Тип корпуса:	F