



Test kit, chloride, model CDS-DT

Товар #: 2580600
RUB Цена: Свяжитесь с Хах
Поставляется в течение 3 дней

The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Спецификации

Smallest increment:	10
Вес с упаковкой:	4
Диапазон измерений:	10 - 10,000 мг/л Cl ⁻
Количество тестов:	100
Наименование метода:	Digital Titrator/ Silver nitrate
Образец тестового набора:	CDS-DT
Платформа:	Digital Titrator
Тип корпуса:	F

Требуемые аксессуары

- Пипетка, переменный объем, 1,0–5,0 мл (Item VBP065)
- Наконечники пипеток для VBP065 (1,0–5,0 мл), 75 шт./уп. (Item VBP068)