



Be Right™



Test kit, water conditioning, model SO-1

Товар #: 225300
RUB Цена: Свяжитесь с Хах
Поставляется в течение 1 недели

Hardness/Soap "Demo" for Sales Aid in Water Conditioning. Use drop-count methods to demonstrate the relationship of hardness to soap consumption. Note: While it provides a quick hardness/soap check, it is less accurate than the hardness titration supplied by Model 15-E (PN 148201) and Model AD-16E (PN 148410) water conditioning test kits.

Uses drop-count methods to demonstrate the relationship of hardness to soap consumption

Provides an easy, economical, and fast hardness/soap check

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Диапазон измерений: 1 drop = 1ppg as CaCO₃
Параметр: Soap Demonstration