

Hardness (Total & Calcium) Test Kit, Model HAC-DT

Товар #:	2063900
RUB Цена:	Связаться с Хах
Доступно	

Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Спецификации

Smallest increment:	0.1-10
Вес с упаковкой:	4
Диапазон измерений:	10 - 4,000 mg/L CaCO ₃
Количество тестов:	100 Ca +100 T
Наименование метода:	Digital Titrator/EDTA
Образец тестового набора:	HAC-DT
Параметр:	Hardness, Total & Calcium - as CaCO ₃
Платформа:	Digital Titrator
Тип корпуса:	F

Что в коробке

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.