



Be Right™



Test kit color disc chlorine dioxide, model CLD-2

Товар #: 2842800
RUB Цена: Свяжитесь с Хах
Доступно

Monitor ClO₂ levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Спецификации

Range:	0 - 2.0 ClO ₂
Smallest increment:	0.1
Вес с упаковкой:	3
Диапазон измерений:	NA
Источник питания:	NA
Количество тестов:	100
Модель:	CLD-2
Наименование метода:	Color Disc/DPD
Образец тестового набора:	CLD-2 ²
Параметр:	Chlorine Dioxide
Платформа:	Color Disc
Тип корпуса:	D