



Be Right™



Test kit, hypochlorite, model CN-HRDT

Товар #: 2687100
RUB Цена: Свяжитесь с Хач
Поставляется в течение 3 дней

Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl₂: 50 - 150 g/L as Cl₂. Contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities

Verify concentration of hypochlorite solutions at time of delivery

Kit contains all required reagents, solutions, and glassware

Спецификации

| | |
|---------------------------|--|
| Smallest increment: | 0.05% |
| Вес с упаковкой: | 3 |
| Диапазон измерений: | 5 - 15 % = 50-150 g/L Cl ₂ |
| Количество тестов: | 100 |
| Наименование метода: | Digital Titrator/ Thiosulfate |
| Номер метода: | 10100 |
| Образец тестового набора: | CN-HRDT |
| Параметр: | Hypochlorite, high range - (Bleach) As Cl ₂ |
| Платформа: | Digital Titrator |
| Тип корпуса: | G |