



Be Right™



Test kit, alkalinity, model AL-TA

Товар #: 2314500
RUB Цена: Свяжитесь с Хах
Доступно

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Спецификации

Smallest increment:	23 gpg
	385
Вес с упаковкой:	1
Диапазон измерений:	23 - 495 gpg
Количество тестов:	100
Наименование метода:	Drop count titration/Sulfuric acid
Образец тестового набора:	AL-TA
Параметр:	Alkalinity - (P) & Total (MO) as CaCO ₃
Платформа:	Drop Count Titrator
Примечание:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Тип корпуса:	B