



Be Right™



## Test kit, alkalinity, model AL-AP

Товар #: 2444300  
RUB Цена: Свяжитесь с Хах  
Доступно

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Спецификации

|                           |   |
|---------------------------|---|
| Smallest increment:       | 0.4 gpg   |
|                           | 1 gpg   |
| Вес с упаковкой:          | 1   |
| Диапазон измерений:       | - 8 gpg   |
| Количество тестов:        | 100   |
| Наименование метода:      | Drop count titration/Sulfuric acid  |
| Образец тестового набора: | AL-AP   |
| Параметр:                 | Alkalinity - (P) & Total (MO) as CaCO <sub>3</sub>  |
| Платформа:                | Drop Count Titrator   |
| Примечание:               | *mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.<br><br>**1 gpg = 17.1 mg/L (MO) = Methyl Orange |
| Тип корпуса:              | B   |