

Fixed volume, Micropipette, 1000 µL

Товар #:	2946010
RUB Цена:	Связаться с Хах
Доступно	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Спецификации

Емкость:	1000 µL
Объем / Размер упаковки:	1 each
Точность, Погрешность:	0,5 %