



Be Right™



## Desiccator, vacuum, 171 mm (6.7") diameter

Товар #: 2238300  
RUB Цена: Свяжитесь с Хах  
Позвоните, чтобы узнать дату доставки

Round polycarbonate top provides lightweight vessel with unmatched transparency. Polypropylene lower section creates support chamber at full vacuum. Stopcock for vacuum applications (1/4-3/8" OD vacuum tubing, sold separately). Holds 28 in Hg (711 mm) for 24 hrs. Silicone O-ring makes grease unnecessary. For room temperature, non-organic solvent, and non-autoclave uses only.

---

### Спецификации

Объем / Размер упаковки: 1 each