



**Be Right™**



## **Flask, Erlenmeyer, glass w/screw cap, 250 mL**

<b>Товар #:</b>	<b>2089646</b>
RUB Цена:	Связаться с Хах
Доступно	

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### **Спецификации**

Емкость:	250 mL
Объем / Размер упаковки:	1 each