



Be Right™



## Beaker, griffin, low form, glass, 250 mL

Товар #:	50046H
RUB Цена:	Связаться с Хах
Доступно	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

---

### Спецификации

Единицы измерения:	each
Емкость:	250 mL
Материал корпуса (датчик):	Glass
Объем / Размер упаковки:	1 each