



Be Right™



Pipet tips, Tensette pipette 1970010, 1.0-10.0 mL, non-sterile, pk/250

Товар #: 2199725

RUB Цена: Свяжитесь с Хах
Доступно

Perform accurate, economical pipetting with Hach's TenSette pipets. Single-channel, air-displacement pipets come in two models: 0.1-1.0 mL volumes (0.1 mL increments) or 1.0-10.0 mL (1.0 mL increments). Volume selection with a simple twist of the turret.

Спецификации

Объем / Размер упаковки: 250 /pk