



Be Right™



## pH paper, 9.0 - 12.0 pH range

Товар #:	38533
RUB Цена:	Связаться с Хах
Доступно	

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

---

### Спецификации

Диапазон измерений:	9.0 - 12.0 pH units
Единицы измерения:	Five 4.6-meter (15') rolls per package