



Be Right™



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 150 mm, circular, 100/pk, qualitative

Товар #:	50658
RUB Цена:	Связаться с Хах
Доступно	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Спецификации

Диаметр:	150 mm
Материал корпуса:	Cellulose
Объем / Размер упаковки:	100 /pk