

Beaker, polypropylene, low form, 250 mL

Товар #:	108046
RUB Цена:	Связаться с Хах
Доступно	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Спецификации

Внешний диаметр:	83.6 mm with spout
Емкость:	250 mL
Материал корпуса:	Polypropylene
Материал корпуса (датчик):	Polypropylene