



Be Right™



Cylinder, graduated, mixing, certified, 100 mL

Товар #: 2636342
RUB Цена: Свяжитесь с Хах
Доступно

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Спецификации

Емкость:	100 mL
Материал корпуса (датчик):	Glass
Метод измерения:	To Contain
Объем / Размер упаковки:	1 each
Сертификаты:	Class A
Тип:	Single Metric, Mixing