

Membrane filters millipore S-pak type HC, 47 mm, 0.7 um pore size, pk/600

Товар #: 2936600
RUB Цена: Свяжитесь с Хах
Доступно

Certified for fecal coliform analysis. Larger pore size (0.7 um) and funnel-like pores creates a less taxing environment for stressed organisms: lessens media loss from evaporation and easier filtration of high particulate samples common to effluents. Whi

Спецификации

Диаметр:	47 mm mm
Материал корпуса:	Mixed-Cellulose Ester
Объем / Размер упаковки:	600
Размер пор:	0.7 μm