



Be Right™



Butterfield's buffered phosphate

Товар #: 2319110
RUB Цена: Свяжитесь с Хах
Позвоните, чтобы узнать дату доставки

Butterfields Buffered Phosphate, 450 mL bottle. pk/10.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Спецификации

Единицы измерения:	10/pk
Объем / Размер упаковки:	450 mL