

Bottle, sampling, presterilized, round, 12/pk

Товар #:	2495012
RUB Цена:	Связаться с Хах
Доступно	

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Works well with P-A Broth Ampules and PathoScreen™ Medium.

Спецификации

Declorinating Agent?:	No
Объем / Размер упаковки:	12/pk